





# PHOTON 4 HEALTH

NEXT LEVEL LIGHT THERAPY


## BODYSLIMMER PRO QUICK START EXPLANATION



You have the choice to treat with manually selected settings or the **SLIM** program (recommended). The above image shows the settings screen. If a code is requested in advance, press 8 repeatedly until the code screen disappears. On the left, at 'Intensity', you can change the light intensity from 1 to 10. We recommend only treating on settings 6, 7, 8, 9 or 10. The thicker the fat layer, the higher the setting. For an average person, 8 is best. Below that, under 'Frequency', you can set the pulsation of the light from 0 to 20 Hz. 'Continuity' means continuous light (0Hz) without pulsation. Finally, under 'Speed', you can adjust the rotation speed of the lasers. We recommend always setting this to 2 or 3. If necessary, adjust the treatment time at 'Remaining' and start the treatment by clicking on the slider next to the SLIM icon 

### The SLIM program

To achieve optimal results, we have developed an automatic 30-minute program that varies and alternates light intensities, rotation speeds and pulsations. In this way, as many fat cells as possible are given the right stimulation for the desired effect. To start this program, press the SLIM button and then the Play icon.

(Pausing and adjusting the lamp height in between treatment remains possible). 

(Optimal treatment distance lasers to body = +-40 cm)

